

Product Datasheet

Human CD200 Receptor 1 (CD200R1) ELISA Kit (orb1203075)

Description	Human CD200 Receptor 1 (CD200R1) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q8TD46
Sensitivity	0.55ng/mL
Entrez	131450
Expiration Date	Please enquire.