

Product Datasheet

Human Chaperonin Containing TCP1, Subunit 2 (CCT2) ELISA Kit (orb1203076)

Description	Human Chaperonin Containing TCP1, Subunit 2 (CCT2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P78371
Sensitivity	0.067ng/mL
Entrez	10576
Expiration Date	Please enquire.