

Product Datasheet

Human Coiled Coil Domain Containing Protein 3 (CCDC3) ELISA Kit (orb1203078)

Description	Human Coiled Coil Domain Containing Protein 3 (CCDC3) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9BQI4
Sensitivity	0.126ng/mL
Entrez	83643
Expiration Date	Please enquire.