

Product Datasheet

Human Bcl2 Like Protein 11 (BCL2L11) ELISA Kit (orb1203084)

Description	Human Bcl2 Like Protein 11 (BCL2L11) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O43521
Sensitivity	0.114ng/mL
Entrez	10018
Expiration Date	Please enquire.