

Product Datasheet

Human ATPase, Na⁺/K⁺ Transporting Beta 4 Polypeptide (ATP1b4) ELISA Kit (orb1203086)

Description	Human ATPase, Na ⁺ /K ⁺ Transporting Beta 4 Polypeptide (ATP1b4) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9UN42
Sensitivity	0.106ng/mL
Entrez	23439
Expiration Date	Please enquire.