

Product Datasheet

Human Aldose Reductase Like Protein 1 (ARL1) ELISA Kit (orb1203091)

Description	Human Aldose Reductase Like Protein 1 (ARL1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	O60218
Sensitivity	0.064ng/mL
Entrez	57016
Expiration Date	Please enquire.