

Product Datasheet

Human Adipocyte Plasma Membrane Associated Protein (APMAP) ELISA Kit (orb1203092)

Description	Human Adipocyte Plasma Membrane Associated Protein (APMAP) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9HDC9
Sensitivity	0.066ng/mL
Entrez	57136
Expiration Date	Please enquire.