

Product Datasheet

Human Aminoacyl tRNA Synthetase Complex Interacting Multifunctional Protein 1 (AIMP1) ELISA Kit (orb1203097)

Description	Human Aminoacyl tRNA Synthetase Complex Interacting Multifunctional Protein 1 (AIMP1) ELISA Kit
Reactivity	Human
Range	125-8000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q12904
Sensitivity	43pg/mL
Entrez	9255
Expiration Date	Please enquire.