

Product Datasheet

Human ATP Binding Cassette Transporter C10 (ABCC10) ELISA Kit (orb1203102)

Description	Human ATP Binding Cassette Transporter C10 (ABCC10) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q5T3U5
Sensitivity	0.061ng/mL
Entrez	89845
Expiration Date	Please enquire.