

Product Datasheet

Rat Tubulin Alpha 4A (TUBa4A) ELISA Kit (orb1203125)

Description	Rat Tubulin Alpha 4A (TUBa4A) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q5XIF6
Sensitivity	0.069ng/mL
Entrez	316531
Expiration Date	Please enquire.