

Product Datasheet

Rat Troponin I Type 1, Slow Skeletal (TNNI1) ELISA Kit (orb1203139)

Description	Rat Troponin I Type 1, Slow Skeletal (TNNI1) ELISA Kit
Reactivity	Rat
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	P13413
Sensitivity	10.7pg/mL
Entrez	29388
Expiration Date	Please enquire.