

Product Datasheet

Rat Solute Carrier Family 30, Member 1 (SLC30A1) ELISA Kit (orb1203167)

Description	Rat Solute Carrier Family 30, Member 1 (SLC30A1) ELISA Kit
Reactivity	Rat
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q62720
Sensitivity	32pg/mL
Entrez	58976
Expiration Date	Please enquire.