

Product Datasheet

Rat Maternal Embryonic Leucine Zipper Kinase (MELK) ELISA Kit (orb1203276)

Description	Rat Maternal Embryonic Leucine Zipper Kinase (MELK) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	D4A6T6
Sensitivity	0.063ng/mL
Entrez	362510
Expiration Date	Please enquire.