

Product Datasheet

Rat Insulin Receptor Substrate 1 (IRS1) ELISA Kit (orb1203298)

Description	Rat Insulin Receptor Substrate 1 (IRS1) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P35570
Sensitivity	0.131ng/mL
Entrez	25467
Expiration Date	Please enquire.