

## Product Datasheet

# Rat Glutamate Receptor, Ionotropic, N-Methyl-D-Aspartate 1 (GRIN1) ELISA Kit (orb1203328)

|                          |                                                                                    |
|--------------------------|------------------------------------------------------------------------------------|
| <b>Description</b>       | Rat Glutamate Receptor,<br>Ionotropic, N-Methyl-D-Aspartate<br>1 (GRIN1) ELISA Kit |
| <b>Reactivity</b>        | Rat                                                                                |
| <b>Range</b>             | 0.156-10ng/mL                                                                      |
| <b>Note</b>              | For research use only                                                              |
| <b>Application notes</b> | sandwich                                                                           |
| <b>Sample Types</b>      | tissue homogenates, cell lysates<br>or other biological fluids.                    |
| <b>Uniprot ID</b>        | <b>P35439</b>                                                                      |
| <b>Sensitivity</b>       | 0.062ng/mL                                                                         |
| <b>Entrez</b>            | <b>24408</b>                                                                       |
| <b>Expiration Date</b>   | Please enquire.                                                                    |