

Product Datasheet

Rat Excitatory Amino Acid Transporter 2 (EAAT2) ELISA Kit (orb1203371)

Description	Rat Excitatory Amino Acid Transporter 2 (EAAT2) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P31596
Sensitivity	0.058ng/mL
Entrez	29482
Expiration Date	Please enquire.