

Product Datasheet

Rat Deiodinase, Iodothyronine, Type I (DIO1) ELISA Kit (orb1203380)

Description	Rat Deiodinase, Iodothyronine, Type I (DIO1) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P24389
Sensitivity	0.059ng/mL
Entrez	25430
Expiration Date	Please enquire.