

Product Datasheet

Rat Cytochrome P450 1A2 (CYP1A2) ELISA Kit (orb1203388)

Description	Rat Cytochrome P450 1A2 (CYP1A2) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P04799
Sensitivity	0.125ng/mL
Entrez	24297
Expiration Date	Please enquire.