

## Product Datasheet

### Rat Carnitine Palmitoyltransferase 1A, Liver (CPT1A) ELISA Kit (orb1203402)

|                          |                                                                      |
|--------------------------|----------------------------------------------------------------------|
| <b>Description</b>       | Rat Carnitine<br>Palmitoyltransferase 1A, Liver<br>(CPT1A) ELISA Kit |
| <b>Reactivity</b>        | Rat                                                                  |
| <b>Range</b>             | 0.312-20ng/mL                                                        |
| <b>Note</b>              | For research use only                                                |
| <b>Application notes</b> | sandwich                                                             |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates or other<br>biological fluids.  |
| <b>Uniprot ID</b>        | <b>P32198</b>                                                        |
| <b>Sensitivity</b>       | 0.107ng/mL                                                           |
| <b>Entrez</b>            | <b>25757</b>                                                         |
| <b>Expiration Date</b>   | Please enquire.                                                      |