

Product Datasheet

Rat Collagen Type IX Alpha 1 (COL9a1) ELISA Kit (orb1203410)

Description	Rat Collagen Type IX Alpha 1 (COL9a1) ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	F1LQ93
Sensitivity	0.52ng/mL
Entrez	305104
Expiration Date	Please enquire.