

Product Datasheet

Rat Collagen Type VI Alpha 1 (COL6a1) ELISA Kit (orb1203411)

Description	Rat Collagen Type VI Alpha 1 (COL6a1) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	D3ZUL3
Sensitivity	0.062ng/mL
Entrez	294337
Expiration Date	Please enquire.