

Product Datasheet

Rat Chymase 1, Mast Cell (CMA1) ELISA Kit (orb1203417)

Description	Rat Chymase 1, Mast Cell (CMA1) ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P50339
Sensitivity	0.57ng/mL
Entrez	25627
Expiration Date	Please enquire.