

Product Datasheet

Rat Ankyrin Repeat Domain Protein 1 (ANKRD1) ELISA Kit (orb1203464)

Description	Rat Ankyrin Repeat Domain Protein 1 (ANKRD1) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q8R560
Sensitivity	0.058ng/mL
Entrez	27064
Expiration Date	Please enquire.