

Product Datasheet

Rat Adrenergic Receptor Beta 1 (ADRB1) ELISA Kit (orb1203475)

Description	Rat Adrenergic Receptor Beta 1 (ADRB1) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P18090
Sensitivity	0.055ng/mL
Entrez	24925
Expiration Date	Please enquire.