

Product Datasheet

Rat Alcohol Dehydrogenase 7 (ADH7) ELISA Kit (orb1203478)

Description	Rat Alcohol Dehydrogenase 7 (ADH7) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P41682
Sensitivity	0.058ng/mL
Entrez	171178
Expiration Date	Please enquire.