

## Product Datasheet

### Porcine Vascular Cell Adhesion Molecule 1 (VCAM1) ELISA Kit (orb1203490)

|                          |                                                               |
|--------------------------|---------------------------------------------------------------|
| <b>Description</b>       | Porcine Vascular Cell Adhesion Molecule 1 (VCAM1) ELISA Kit   |
| <b>Reactivity</b>        | Porcine                                                       |
| <b>Range</b>             | 15.6-1000ng/mL                                                |
| <b>Note</b>              | For research use only                                         |
| <b>Application notes</b> | sandwich                                                      |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates or other biological fluids. |
| <b>Uniprot ID</b>        | <b>Q28939</b>                                                 |
| <b>Sensitivity</b>       | 6.7ng/mL                                                      |
| <b>Entrez</b>            | <b>396925</b>                                                 |
| <b>Expiration Date</b>   | Please enquire.                                               |