

Product Datasheet

Bovine Transthyretin (TTR) ELISA Kit (orb1203497)

Description	Bovine Transthyretin (TTR) ELISA Kit
Reactivity	Bovine
Range	7.81-500ng/mL
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O46375
Sensitivity	2.8ng/mL
Entrez	280948
Expiration Date	Please enquire.