

Product Datasheet

Bovine TEK Tyrosine Kinase, Endothelial (Tie2) ELISA Kit (orb1203514)

Description	Bovine TEK Tyrosine Kinase, Endothelial (Tie2) ELISA Kit
Reactivity	Bovine
Range	6.25-400ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q06807
Sensitivity	1.56ng/mL
Entrez	280939
Expiration Date	Please enquire.