

Product Datasheet

Fish Tyrosine Kinase With Immunoglobulin Like And EGF Like Domains Protein 1 (Tie1) ELISA Kit (orb1203515)

Description	Fish Tyrosine Kinase With Immunoglobulin Like And EGF Like Domains Protein 1 (Tie1) ELISA Kit
Reactivity	Zebrafish
Range	7.81-500ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	073792
Sensitivity	3.2ng/mL
Entrez	30746
Expiration Date	Please enquire.