

## Product Datasheet

# Bovine Peroxisome Proliferator Activated Receptor Gamma (PPAR $\gamma$ ) ELISA Kit (orb1203544)

|                          |                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Bovine Peroxisome Proliferator<br>Activated Receptor Gamma<br>(PPAR $\gamma$ ) ELISA Kit                      |
| <b>Reactivity</b>        | Bovine                                                                                                        |
| <b>Range</b>             | 0.312-20ng/mL                                                                                                 |
| <b>Note</b>              | For research use only                                                                                         |
| <b>Application notes</b> | sandwich                                                                                                      |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates, cell lysates, cell<br>culture supernates or other<br>biological fluids. |
| <b>Uniprot ID</b>        | <b>O18971</b>                                                                                                 |
| <b>Sensitivity</b>       | 0.115ng/mL                                                                                                    |
| <b>Entrez</b>            | <b>281993</b>                                                                                                 |
| <b>Expiration Date</b>   | Please enquire.                                                                                               |