

Product Datasheet

Bovine Oncostatin M (OSM) ELISA Kit (orb1203552)

Description	Bovine Oncostatin M (OSM) ELISA Kit
Reactivity	Bovine
Range	62.5-4000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P53346
Sensitivity	15.8pg/mL
Entrez	319086
Expiration Date	Please enquire.