

Product Datasheet

Porcine Matrix Metalloproteinase 3 (MMP3) ELISA Kit (orb1203568)

Description	Porcine Matrix Metalloproteinase 3 (MMP3) ELISA Kit
Reactivity	Porcine
Range	4.68-300ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	C9VZX5
Sensitivity	1.8ng/mL
Entrez	396769
Expiration Date	Please enquire.