

Product Datasheet

Bovine Lactoferrin (LTF) ELISA Kit (orb1203585)

Description	Bovine Lactoferrin (LTF) ELISA Kit
Reactivity	Bovine
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, saliva, milk, cell culture supernates or other biological fluids.
Uniprot ID	P24627
Sensitivity	0.71ng/mL
Entrez	280846
Expiration Date	Please enquire.