

Product Datasheet

Porcine Leukemia Inhibitory Factor (LIF) ELISA Kit (orb1203587)

Description	Porcine Leukemia Inhibitory Factor (LIF) ELISA Kit
Reactivity	Porcine
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9GKZ8
Sensitivity	0.047ng/mL
Entrez	399503
Expiration Date	Please enquire.