

Product Datasheet

Bovine Leukemia Inhibitory Factor (LIF) ELISA Kit (orb1203588)

Description	Bovine Leukemia Inhibitory Factor (LIF) ELISA Kit
Reactivity	Bovine
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q27956
Sensitivity	11.5pg/mL
Entrez	280840
Expiration Date	Please enquire.