

Product Datasheet

Bovine Galectin 4 (GAL4) ELISA Kit (orb1203623)

Description	Bovine Galectin 4 (GAL4) ELISA Kit
Reactivity	Bovine
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q3T0D6
Sensitivity	0.112ng/mL
Entrez	614804
Expiration Date	Please enquire.