

Product Datasheet

Porcine Factor Related Apoptosis Ligand (FASL) ELISA Kit (orb1203629)

Description	Porcine Factor Related Apoptosis Ligand (FASL) ELISA Kit
Reactivity	Porcine
Range	23.4-1500pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9BEA8
Sensitivity	5.2pg/mL
Entrez	396726
Expiration Date	Please enquire.