

Product Datasheet

Bovine C Reactive Protein (CRP) ELISA Kit (orb1203642)

Description	Bovine C Reactive Protein (CRP) ELISA Kit
Reactivity	Bovine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	C4T8B4
Sensitivity	6.2pg/mL
Entrez	527553
Expiration Date	Please enquire.