

Product Datasheet

Bovine Creatine Kinase, Muscle (CKM) ELISA Kit (orb1203646)

Description	Bovine Creatine Kinase, Muscle (CKM) ELISA Kit
Reactivity	Bovine
Range	12.5-800U/L
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9XSC6
Sensitivity	4.65U/L
Entrez	286822
Expiration Date	Please enquire.