

Product Datasheet

Bovine Cadherin, Epithelial (CDHE) ELISA Kit (orb1203648)

Description	Bovine Cadherin, Epithelial (CDHE) ELISA Kit
Reactivity	Bovine
Range	0.937-60ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q6R8F2
Sensitivity	0.23ng/mL
Entrez	282637
Expiration Date	Please enquire.