

Product Datasheet

Bovine Cluster Of Differentiation 40 Ligand (CD40L) ELISA Kit (orb1203650)

Description	Bovine Cluster Of Differentiation 40 Ligand (CD40L) ELISA Kit
Reactivity	Bovine
Range	0.625-40ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P51749
Sensitivity	0.18ng/mL
Entrez	282387
Expiration Date	Please enquire.