

Product Datasheet

Bovine Apolipoprotein H (APOH) ELISA Kit (orb1203664)

Description	Bovine Apolipoprotein H (APOH) ELISA Kit
Reactivity	Bovine
Range	3.12-200ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P17690
Sensitivity	1.87ng/mL
Entrez	281006
Expiration Date	Please enquire.