

## Product Datasheet

### Bovine Angiopoietin 1 (ANGPT1) ELISA Kit (orb1203674)

|                          |                                                                     |
|--------------------------|---------------------------------------------------------------------|
| <b>Description</b>       | Bovine Angiopoietin 1 (ANGPT1)<br>ELISA Kit                         |
| <b>Reactivity</b>        | Bovine                                                              |
| <b>Range</b>             | 0.625-40ng/mL                                                       |
| <b>Note</b>              | For research use only                                               |
| <b>Application notes</b> | sandwich                                                            |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates or other biological<br>fluids. |
| <b>Uniprot ID</b>        | <b>O18920</b>                                                       |
| <b>Sensitivity</b>       | 0.195ng/mL                                                          |
| <b>Entrez</b>            | <b>282140</b>                                                       |
| <b>Expiration Date</b>   | Please enquire.                                                     |