

Product Datasheet

Bovine Anti-Mullerian Hormone (AMH) ELISA Kit (orb1203676)

Description	Bovine Anti-Mullerian Hormone (AMH) ELISA Kit
Reactivity	Bovine
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P03972
Sensitivity	0.117ng/mL
Entrez	280718
Expiration Date	Please enquire.