

Product Datasheet

Porcine Angiotensin I Converting Enzyme (ACE) ELISA Kit (orb1203685)

Description	Porcine Angiotensin I Converting Enzyme (ACE) ELISA Kit
Reactivity	Porcine
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	A1YQ85
Sensitivity	10.7pg/mL
Entrez	613133
Expiration Date	Please enquire.