

Product Datasheet

Bovine Angiotensin I Converting Enzyme (ACE) ELISA Kit (orb1203686)

Description	Bovine Angiotensin I Converting Enzyme (ACE) ELISA Kit
Reactivity	Bovine
Range	6.25-400ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q9BDG1
Sensitivity	2.15ng/mL
Entrez	509484
Expiration Date	Please enquire.