

Product Datasheet

Porcine Fibroblast Growth Factor 21 (FGF21) ELISA Kit (orb1203695)

Description	Porcine Fibroblast Growth Factor 21 (FGF21) ELISA Kit
Reactivity	Porcine
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	C6JWG4
Sensitivity	13.3pg/mL
Entrez	100302504
Expiration Date	Please enquire.