

Product Datasheet

Rat p53 Upregulated Modulator Of Apoptosis (PUMA) ELISA Kit (orb1203734)

Description	Rat p53 Upregulated Modulator Of Apoptosis (PUMA) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q80ZG6
Sensitivity	0.047ng/mL
Entrez	317673
Expiration Date	Please enquire.