

Product Datasheet

Rat Protein Kinase, AMP Activated Beta 1 (PRKAb1) ELISA Kit (orb1203737)

Description	Rat Protein Kinase, AMP Activated Beta 1 (PRKAb1) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P80386
Sensitivity	0.060ng/mL
Entrez	83803
Expiration Date	Please enquire.