

## Product Datasheet

### Rat Lipase, Pancreatic (PL) ELISA Kit (orb1203742)

|                          |                                              |
|--------------------------|----------------------------------------------|
| <b>Description</b>       | Rat Lipase, Pancreatic (PL)<br>ELISA Kit     |
| <b>Reactivity</b>        | Rat                                          |
| <b>Range</b>             | 78.1-5000pg/mL                               |
| <b>Note</b>              | For research use only                        |
| <b>Application notes</b> | sandwich                                     |
| <b>Sample Types</b>      | serum, plasma or other<br>biological fluids. |
| <b>Uniprot ID</b>        | <b>P27657</b>                                |
| <b>Sensitivity</b>       | 37pg/mL                                      |
| <b>Entrez</b>            | <b>25702</b>                                 |
| <b>Expiration Date</b>   | Please enquire.                              |